



## Polytetrafluoroethylene Felt Fabric

Offered By TDA

For Medical & Food Applications

Smooth Texture  
Medium & High Densities  
Fibrous  
Conformable

Consider PTFE Felts for...

- Pledgets
- Buttresses
- Plugs
- Washers
- Seals

Standard Sizes Available  
4" X 4" (10 cm X 10 cm)  
6" X 6" (15 cm X 15 cm)

Custom sizes available upon request.

Properties	PTFENF900S	PTFENF900ST	PTFENF1800S	PTFENF1800ST
Weight, g/m <sup>2</sup> (ASTM D-3776)	900	900	1800	1800
Thickness, mm (ASTM D-5729)	1.50	1.00	1.85	1.50
Strength MD, N/2.5 cm (ASTM D-5034)	601	601	1442	1442
Strength CMD, N/2.5 cm (ASTM D-5034)	890	890	1146	1146
Mullen Burst, kPa (ASTM D-751)	2310	2310	4833	4833
Apparent Density, g/cc (Calculated)	0.60	0.90	0.97	1.20
Porosity, % (Void Volume)	73	59	56	45
Relative Density	Medium	High	High	High

**DISCLAIMER:** The information supplied in this document is for guidance only and should not be construed as a warranty. All implied warranties are expressly disclaimed, including without limitation any warranty of merchantability and fitness for use. All users of the material are responsible for assuring that it is suitable for their needs, environment and end use. All data is subject to change as TDA deems appropriate.